REF 925703

05.22

de

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05.22

en

NANOCOLOR® Analytical Quality Control

NANOCONTROL NANOCHECK 2.0

Method:

Test solutions for the determination of photometric accuracy and linearity (secondary standards in accordance to ISO 9001)

Suitable for *NANOCOLOR®* photometers (400 D / 500 D / UV_{VIS} / VIS / $VIS II$ / $UV_{VIS II}$ / Advance) and PF-12 / PF-12^{Plus} / PF-3

Factor:	1.000
Wavelength (HW = 5 – 12 nm):	345 – 800 nm
Reaction temperature:	20 – 25 °C

Contents of reagent set:

- 2 tubes with test solution 1.1
- 2 tubes with test solution 1.2
- 2 tubes with test solution 1.3
- 2 tubes with test solution 2.1
- 2 tubes with test solution 2.2
- 2 tubes with test solution 2.3
- 1 tube with blank value "NULL"
- 1 instructions for use
- 1 lot-specific certificate

Note:

The test solutions have been checked in a reference photometer monitored with primary standards (NIST-Standards). The results are documented.

Hazard warning:

Information regarding safety can be found on the box' label and in the safety data sheet. You can download the SDS from www.mn-net.com/SDS.

Procedure:

1. Call up the extinction program (e.g. method 903 on photometer *NANOCOLOR®* 500 D)
2. Choose the wavelength you want to check from the evaluation sheet*
3. Measure the blank value "NULL"
4. Insert different test solutions one after another and measure
5. Record the result in the evaluation sheet

Note: * Device-specific evaluation sheets can be downloaded on www.mn-net.com/NANOCHECK. Tubes must be clean. If necessary, clean them with a dry and clean cloth. Do not use scratched tubes.

Please note when using *NANOCOLOR®* PC software for spectrophotometer:

Test solution 1 $\triangleq$ Test solution 1.1; Test solution 2 $\triangleq$ Test solution 1.3

Interpretation:

The measuring values of the different wavelengths with all test solutions are recorded on the evaluation sheet and have to be compared with the given nominal value.

If the measuring values are **within the confidence interval**, the wavelength settings and measuring geometry of the photometer are in order.

If the measuring values are **out of the confidence interval**, please proceed as follows:

1. Check the power supply
2. Change the test solutions (the reagent set contains every test solutions twice)
3. Check expiry date of the reagent set; if necessary, use a new pack
4. Check cleanliness of the tubes
5. Check if tubes are scratched
6. Calibrate the photometer and repeat the NANOCHECK procedure
7. If necessary, contact your distributor or MACHEREY-NAGEL

Disposal:

Information regarding disposal can be found in the safety data sheet. You can download the SDS from www.mn-net.com/SDS.

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