REF 918979

06.20

de

NANOCOLOR®

NanOx N

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NANOCOLOR®

NanOx N

Solid reagents for oxidative decomposition

to determine total nitrogen TN_b

Contents:

- 30 g *NanOx* N decomposition reagent
- 1 spoon 85 mm orange
- 25 g *NanOx* N compensation reagent
- 1 spoon 85 mm black

Hazard warning:

Information on hazards can be found on the outer label and on the safety data sheet. The safety data sheet can be downloaded from www.mn-net.com/SDS.

Method:

Oxidation of all inorganic and organic nitrogen-containing substances to form nitrate. Residues of peroxide and eventually oxidized chromium(VI), which remain after oxidation and interfere with the determination of nitrate, are eliminated by a compensation reagent.

Field of application:

Decomposition in a heating block

The decomposition in a heating block at 120 °C features a lower oxidation potential than the microwave decomposition. For mainly municipal wastewaters this method can be used, if the matrix is constant over a longer time period. In certain intervals the applicability of the method should be checked against a reference procedure (e.g. microwave decomposition). For industrial wastewaters, however, nitrogen compounds which are difficult to oxidize, frequently are not or only partially decomposed in a heating block.

Pressurized microwave decomposition

This method offers the advantage of considerable time-saving and convenient operation compared to wet decompositions in the heating block. The recovery of organic components, which are difficult to decompose, is very high.

Determination of total nitrogen TN_b (total N):

Ranges:	0.5–22.0 mg/L N	method 0831
	0.3–8.0 mg/L N	method 0651
	3–60 mg/L N	method 0921
	5–220 mg/L N	method 0881

The inorganic components ammonium, nitrite and nitrate are determined together with the organic nitrogen compounds such as aminoacids, urea, complexing agents etc. The determination results as nitrate nitrogen with test 0-64 *NANOCOLOR*® Nitrate 50.

Remarks:

The pH value of the sample to be decomposed must be between pH 5 and 9, if necessary adjust with sodium hydroxide solution or sulphuric acid. Nitrogen concentrations above the double measuring range can simulate results within the measuring range and thus cause a wrong evaluation. Dilute the sample until the measured value is within the measuring range. For waters of unknown concentrations we recommend that you perform the test with very different dilutions (e.g. 1+9, 1+99) until the last dilution confirms the previous value. For samples which consume large amounts of oxidizing substances (e.g. for COD values above 1000 mg/L O₂), decomposition can be incomplete, or the safety diaphragm in the pressure vessel can be destroyed. In such cases repeat the decomposition with a diluted sample solution.

Determination of total phosphorus (total P) and total metal:

For the determination of total phosphorus (sum of all phosphorus compounds such as ortho-phosphates, polyphosphates and organophosphates) as well as for metals bonded as complexes or for metals, which are present in an oxidation state which is not covered by an analytical procedure (e.g. chromium(III)), you have to perform the decomposition with reagent set *NANOCOLOR*® *NanOx* Metall (REF 918978).

Procedure for decomposition in the heating block:

Requisite accessories:

NANOCOLOR® heating block
NANOCOLOR® 16 mm OD test tubes (REF 91680)

Pipet 5.0 mL sample solution into an **empty** reaction tube 16 mm OD, add 1 level orange measuring spoon *NanOx* N decomposition reagent, close and shake thoroughly. Place the reaction tube into the heating block and heat 30 min at 120 °C or 1 h at 100 °C. Remove tube from the heating block, shake lightly and let cool. Open the reaction tube and add 1 level black measuring spoon *NanOx* N compensation reagent, close and shake thoroughly. The contents can now be used as sample solution for determination of total nitrogen TN_b with test 0-64 *NANOCOLOR*® Nitrate 50.

Alternative procedure for microwave decomposition:

Requisite accessories:

Microwave oven with digital display (time setting in seconds), PTFE pressure vessel 20 mL for decomposition and safety expansion reservoir (REF 91640)

Pipet 10 mL sample solution into the decomposition vessel, add 2 level orange measuring spoons *NanOx* N decomposition reagent, close and shake thoroughly. Place pressure vessel on the outer rim of the turntable in the microwave oven and heat 23 s at 900 VA or 28 s at 750 VA (always use the highest setting of the respective instrument). Remove vessel from the microwave oven and let cool for about 10 min. Turn the pressure vessel on its head once and then open it with caution. Add 2 level black measuring spoons *NanOx* N compensation reagent, close and shake thoroughly. The contents can now be used as sample solution for determination of total nitrogen TN_b with test 0-64 *NANOCOLOR*® Nitrate 50.

Dilutions:

Sample 1+4 diluted	Range 3–60 mg/L N
Heating block:	Pipet 1.0 mL test sample and 4.0 mL dist. water into a 14 mm ID reaction tube
Microwave:	Pipet 2.0 mL test sample and 8.0 mL dist. water into the decomposition vessel
Sample 1+9 diluted	Range 5–220 mg/L N
Heating block:	Pipet 0.5 mL test sample and 4.5 mL dist. water into a 14 mm ID reaction tube
Microwave:	Pipet 1.0 mL test sample and 9.0 mL dist. water into the decomposition vessel

Interferences:

- a) **The safety diaphragm ruptures**
- 1- Material fatigue
→ replace safety diaphragm and sealing disc
 - 2- Sample develops gases → dilute sample
 - 3- Delays in boiling
→ use special *NanOx* N boiling chips
 - 4- Wrong time setting → adjust time
- b) **Discoloration of the Nitrate 50 test tube (can simulate positive results)**
- 1- Decomposition incomplete
Note: slightly nitrous odour when opening the decomposition vessel or decomposition solution not clear
→ dilute into the measuring range, if necessary decompose twice
 - 2- other interferences, which cannot be eliminated by decomposition with *NanOx* N

Time adjustment for different microwave ovens:

For decomposition with *NanOx* N the microwave times have to be adjusted for different instruments. All instruments should be operated at the highest setting, and the time is adjusted to the second.

Example for calculation:

Microwave energy 900 VA → 24 s heating time

Microwave energy 750 VA → X s heating time

X = 900/750 x 24 s = 28.8 s → 29 s

Please note that the actual power of the microwave oven may deviate from the nominal value. In these cases test and adjust the time as described below.

Correction and optimization of microwave times:

The output of microwave ovens can fluctuate even between instruments of the same model or it can change with time. For the decomposition with *NanOx* N, however, an exact dosing of energy is required. In order to optimize the microwave time, always choose the highest setting of the respective instrument (max. 1000 VA), because this is not interrupted by time intervals. For a microwave oven with 900 VA output the microwave time will be about 24 s. The following procedure is a very reliable and simple method to determine the microwave energy:

Fill the included Erlenmeyer flask 100 mL with exactly 100 mL water (of ambient temperature) and measure the temperature with a thermometer (e.g. 21 °C). Now place the filled Erlenmeyer flask (without thermometer) on the outer rim of the microwave turntable (same position as the decomposition vessel) and heat for the calculated time. Immediately remove the Erlenmeyer flask and measure the water temperature after slight stirring (e.g. 53 °C). The temperature difference must

be 30 °C. In our example the difference is 53 - 21 = 32 °C. The heating time has to be corrected for this deviation:

Nominal difference/actual difference x time = corrected time
30 / 32 x 24 s = 22.5 s → 23 s

Check this time setting again following the above procedure.

Note: For microwave decomposition of total nitrogen with *NanOx* N you may lower the heating time by about 5 % to prevent destruction of the safety diaphragm (e.g. 22 s in the above example).

Control of decomposition parameters:

With the *NANOCONTROL* multistandards you can easily check the microwave decomposition using the determined time setting. The multistandards contain bonded nitrogen, ammonium and nitrate. Analyze for these parameters prior to decomposition.

General remarks:

Always place the decomposition vessel on the same position of the microwave turntable. For easier handling we recommend that you fill the safety expansion reservoir with filter flocs or pulp, covered with a filter circle, and remove the exhaust tubing from the safety expansion reservoir. After 5 decompositions the filter flocs or pulp have to be replaced. If the safety diaphragm ruptures, always dry the safety reservoir prior to the next run, or replace the filter flocs or pulp. For sample solutions which tend to delays in boiling, we recommend using PTFE boiling chips. After decomposition remove the decomposition solution from the vessel. When the vessel is heavily soiled, or after continued use, clean the decomposition vessel with a bottle brush and dist. water. We recommend that from time to time you perform a blank decomposition with dist. water and *NanOx* N decomposition reagent, but without addition of compensation reagent. As an alternative, you may fill the vessel with dist. water and some drops of hydrogen peroxide 30 % and let it stand for some hours.

REF	
918979	<i>NANOCOLOR® NanOx</i> N solid reagents 30 g decomposition reagent plus measuring spoon 25 g compensation reagent plus measuring spoon
91640	PTFE pressure vessel 20 mL for decomposition with safety expansion reservoir, set of sealing discs and safety diaphragms
91643 91640.1 91646 91644 91645 91628 91680 481100 481110	PTFE pressure vessel 20 mL for decomposition PFA safety expansion reservoir PFA tubing with fitting pack of 20 PFA sealing discs pack of 10 PTFE safety diaphragms pack of 3 PTFE boiling chips pack of 20 reaction tubes 16 mm OD filter flocs, pack of 500 g filter flocs, pack of 1000 g

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MACHEREY-NAGEL GmbH & Co. KG
France : MACHEREY-NAGEL SAS
MACHEREY-NAGEL SAS