REF 91808

04.23

de

NANOCOLOR® Aufschluss-Set**Methode:**

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en

NANOCOLOR® Crack Set**Method:**

Dissolution and decomplexation of heavy metals with sulfuric acid and potassium peroxodisulfate

Contents of reagent set:

100 mL Crack Set R1
17 g Crack Set R2
100 mL Crack Set R3
1 measuring spoon 85 mm

Hazard warning:

Information regarding safety can be found on the box' label and in the safety data sheet. You can download the SDS from www.mn-net.com/SDS.

Procedure:

Requisite accessories: heating block *NANOCOLOR® VARIO C2 M* (REF 919350.1), decomposition apparatus (REF 91629), piston pipette with tips

Switch on heating block, set to 100 °C and 1:00 h.

Fill decomposition tube with
10 mL homogenized sample, add
1 mL R1 and

1 level measuring spoon R2, shake slightly, attach absorption tube and place decomposition tube into the heating block, press **START** key.

After 1 h remove decomposition tube from heating block, cool, add
1 mL R3 and mix. *The pH value should be between a pH of 2 and 5, otherwise add more or less R3.*

A. Tube tests

Name	Test	Method
Lead 5	0-09	0093
Cadmium 2	0-14	0143
Iron 3	0-37	0373
Zinc 6	0-42	0423
Copper 5	0-53	0533

Name	Test	Method
Copper 7	0-54	0543
Nickel 7	0-61	0613
Nickel 4	0-71	0713
Zinc 4	0-96	0963

Measurement:

For MACHEREY-NAGEL photometers see manual.

Photometers of other manufacturers:

Using test tubes please notice the dilution of the decomposed solution! **Multiply results by 1.2!**

B. Standard tests

Name	Test
Iron*	1-36
Cobalt	1-51
Copper	1-53

Name	Test
Nickel	1-62
Zink	1-95

* An exact determination of iron requires a reagent blank value, because even analytical grade chemicals may contain traces of iron.

Measurement:

For MACHEREY-NAGEL photometers see manual. **Multiply results by 2!**

Photometers of other manufacturers:

Using standard tests please notice the dilution of the decomposed solution! **Multiply results by 2!**

C. Extraction methods

Pour decomposed solution into a separation funnel 100 mL and rinse decomposition tube with about 40 mL distilled water (corresponds to the 50 mL sample solution mentioned in the procedures).

Suited for: Test 1-10 Lead, Test 1-13 Cadmium.

Measurement:

For MACHEREY-NAGEL photometers see manual. **Multiply results by 5!**

Photometers of other manufacturers:

Using extraction methods please notice the dilution of the decomposed solution! **Multiply results by 5!**

Disposal:

Information regarding disposal can be found in the safety data sheet. You can download the SDS from www.mn-net.com/SDS.

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