


# visicolor® ECO Sulfid


## Sulphide / Sulfure / Sulfuro

Method(e) / Método


5941 0.05–0.80 mg/L S<sup>2-</sup>PF-12  
620 nmPF-12<sup>Plus</sup>  
620 nmAdvance  
670 nm

Null mit Probe  
Blank with sample  
Blanc avec échantillon  
En blanco con muestra


Geeignet für Meerwasser  
Suitable for seawater  
Convient à l'eau de mer  
Adecuado para el agua de mar




Küvetten spülen  
Rinse test tube  
Rincer les éprouvettes  
Lavar los tubos




Probe  
Sample  
Echantillon  
Muestra




1 S-1




Umschwenken  
Shake gently  
Secouer légèrement  
Mezclar volteándolo




Warten  
Wait  
Attendre  
Espere




5 S-2



5 S-3




Schütteln  
Shake  
Agiter  
Agitar



Säubern  
Clean  
Nettoyer  
Limpiar



Warten  
Wait  
Attendre  
Espere



Messung  
Measurement  
Mesure  
Medición