

MACHEREY-NAGEL

Column racks for anion exchange chromatography



NucleoBond® Racks for plasmid purification

- Racks for Mini to Giga preps
- New height adjustable NucleoBond® PC 10000 Rack
- Convenient plasmid purification

NucleoBond® Racks for plasmid purification

NucleoBond® Racks

One size does not fit all – Get the right rack for your NucleoBond® column!

Mini to Giga – Get the right rack for convenient gravity flow processing and elution.



Follow the link to our new workflow video and learn how simple it is to process our NucleoBond® PC 10000 and 10000 EF Giga Kits in combination with our new height adjustable NucleoBond® PC 10000 Rack.



Ordering information

Product	Specifications	REF	Format
 NucleoBond® Rack Small	For use with NucleoBond® AX 20 Columns, RNA Soil Mini Columns 	740562	Mini 
 NucleoBond® Rack Large	For use with NucleoBond® AX 100, AX 500, AX 2000, BAC 100, NucleoBond® HMW and Xtra Midi Columns, plastic washers are needed for all columns 	740563	Midi Maxi  
 NucleoBond® Xtra Combi Rack	For use with NucleoBond® Xtra Midi, Xtra Maxi, AX 500, Xtra BAC, BAC 100, NucleoBond® HMW and NucleoBond® RNA Columns 	740415	Midi Maxi  
 NucleoBond® Smart Rack	For use with NucleoBond® Xtra Midi, Xtra Maxi, AX 500, AX 2000, BAC 100, NucleoBond® HMW and NucleoBond® RNA Columns 	740413	Midi Maxi  
 NucleoBond® PC 10000 Rack	For use with NucleoBond® AX 10000 Columns 	740599	Giga 

www.mn-net.com

MACHERY-NAGEL



MACHERY-NAGEL GmbH & Co. KG
Valencienner Str. 11
52355 Düren · Germany

DE Tel.: +49 24 21 969-0 info@mn-net.com
CH Tel.: +41 62 388 55 00 sales-ch@mn-net.com
FR Tel.: +33 388 68 22 68 sales-fr@mn-net.com
US Tel.: +1 888 321 62 24 sales-us@mn-net.com

