


# visicolor® ECO Cyanid


## Cyanide / Cyanure / Cianuro

Method(e) / Método


5221 0.01–0.20 mg/L CN<sup>-</sup>PF-12  
585 nmPF-12<sup>Plus</sup>  
585 nmAdvance  
605 nm

Null mit Probe  
Blank with sample  
Blanc avec échantillon  
En blanco con muestra


Geeignet für Meerwasser nach Verdünnung 1+3  
Suitable for sea water after dilution 1+3  
Convient à l'eau de mer après dilution 1+3  
Adecuado para el agua de mar tras su dilución 1+3




Küvetten spülen  
Rinse test tube  
Rincer les éprouvettes  
Lavar los tubos




Probe  
Sample  
Echantillon  
Muestra




5 Δ CN-1




Schütteln  
Shake  
Agiter  
Agitar




1 — CN-2




Umschwenken  
Shake gently  
Secouer légèrement  
Mezclar volteándolo



5 Δ CN-3




Schütteln  
Shake  
Agiter  
Agitar



Säubern  
Clean  
Nettoyer  
Limpiar



Warten  
Wait  
Attendre  
Espere



Messung  
Measurement  
Mesure  
Medición