

Filtration

Filter papers for sugar refineries



MACHEREY-NAGEL



Filter papers for sugar refineries

Control of incoming sugar beets Reference



Determination of sugar content
automized filtration in beetanalyzing systems

MN 672 reels, 1000 m x 240 mm

100 672 001

Determination of sugar content, alternatively to reels



MN 614 filter circles, 240 mm Ø
MN 614 filter circles ¼ folded, 220 mm
MN 620 filter sheets, 580 x 580 mm

427 024

400 614 022 G

141 010

Quality control



Determination of sugar content and colour

MN 614 filter circles, 240 mm Ø
MN 614 filter circles ¼ folded, 220 mm

427 024

400 614 022 G

Weighing of thick liquids



Parchment sheets MN 40/25
100 x 100 mm, pack of 1000 sheets

194 000

Filtration of suspended solids
from sugar solutions



Membrane filters PORAFIL® CA, 0,45 µm, 50 mm Ø

680000 450 50

Sewage plants



Determination of dry residue / ash content
for suction filter 90 mm

MN 640 w filter circles, 90 mm Ø
MN 640 d filter circles, 90 mm Ø

202 009

205 009



Determination of dry residue / ash content
for suction filter 110 mm

MN 640 w filter circles, 110 mm Ø
MN 640 d filter circles, 110 mm Ø

202 011

205 011

Your local distributor:



With kind permission from Strube-Dieckmann, 31688 Nienstädt, Germany

KATEN100067Zuckerfabriken em1/20/7.2006 PD
Printed in Germany

Further information and products: www.mn-net.com

MACHEREY-NAGEL

MACHEREY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6-8 · 52355 Düren · Germany

Germany and international:

Tel.: +49 (0) 24 21 96 90

Fax: +49 (0) 24 21 96 91 99

e-mail: sales-de@mn-net.com

Switzerland:

Tel.: +41 (0) 62 388 55 00

Fax: +41 (0) 62 388 55 05

e-mail: sales-ch@mn-net.com

France:

Tel.: +33 (0) 3 88 68 22 68

Fax: +33 (0) 3 88 51 76 88

e-mail: sales-fr@mn-net.com

USA:

Tel.: +1 610 559 98 45

Fax: +1 610 559 98 78

e-mail: sales-us@mn-net.com

