

NucleoMag® SEP Maxi



Dear valued customer,

Thank you for choosing MACHEREY-NAGEL for your application.

This product is for research use only.

Important information regarding product components, specifications, safety instructions, and processing protocols, can all be found on the product website and accessed easily via the QR code.

[QR-Code product website](#)



[qr.mn-net.com/qr/\(241\)744902](https://qr.mn-net.com/qr/(241)744902)

Contents

Product	REF	Pack of
NucleoMag® SEP Maxi	744902	1

Product description

The NucleoMag® SEP Maxi magnetic separator is designed for the rapid manual purification of nucleic acids using the NucleoMag® kits. The NucleoMag® kits are based on MACHEREY-NAGEL's magneticbead technology and allow rapid simultaneous purification of highquality nucleic acids.

The NucleoMag® SEP Maxi is used with 50 mL tubes and carries 4 tube positions.

Specifications

The NucleoMag® SEP Maxi consists of a Polypropylen (PP) base. Do not autoclave. If necessary, clean the NucleoMag® SEP Maxi by rinsing it with water or by wiping off dirt and liquids with a towel. NucleoMag® SEP Maxi is resistant to diluted acetic acid, ethanol, guanidine-HCl, NaCl, NaOH, and SDS. However, prolonged exposure of the separator to these chemicals should be avoided.

Ordering information

Product	REF	Pack of
NucleoMag® SEP Magnetic separator for 96-well plates	744900	1
NucleoMag® SEP Mini Magnetic separator for 1.5 mL or 2 mL tubes	744901	1
NucleoMag® SEP Maxi Magnetic separator for 50 mL tubes	744902	1

Safety instructions

When working with this product, wear suitable protective clothing (e.g., lab coat, disposable gloves, and protective goggles).

[QR-Code Safety Data Sheet](#)

Please find a digital version of the safety data sheets by following the link below:



www.mn-net.com/sds

Product use restriction / warranty

All MACHEREY-NAGEL products are designed for their intended use only. They are not intended to be used for any other purpose. The description of the intended use of the products can be found in the original MACHEREY-NAGEL product leaflets. Before using our products, please observe the instructions for use and the safety instructions from the respective Material Safety Data Sheet of the product.

This MACHEREY-NAGEL product is carrying documentation stating specifications and other technical information. MACHEREY-NAGEL warrants to meet the stated specifications. The provided warranty is limited to the data specifications and descriptions as given in the original MACHEREY-NAGEL literature. No other statements or representations, written or oral, by MACHEREY-NAGEL's employees, agents or representatives, except written statements signed by a duly authorized officer of MACHEREY-NAGEL are authorized. They should not be relied upon by the customer and are not a part of a contract of sale or of this warranty.

Liability for all possible damages that occur in any connection with our products is limited to the utmost minimum as stated in the general business terms and conditions of MACHEREY-NAGEL in their latest edition which can be taken from the company's website. MACHEREY-NAGEL does not assume any further warranty.

Products and their application are subject to change. Therefore, please contact our Technical Service Team for the latest information on MACHEREY-NAGEL products. You may also contact your local distributor for general scientific information. Descriptions in MACHEREY-NAGEL literature are provided for informational purposes only.

Last updated: 08/2022, Rev. 04

Please contact:

MACHEREY-NAGEL GmbH & Co. KG

Tel.: +49 24 21 969-333

support@mn-net.com

Contact MN

Germany and international

MACHEREY-NAGEL GmbH & Co. KG

Valenciennner Str. 11 · 52355 Düren · Germany

Tel.: +49 24 21 969-0

Toll-free: 0800 26 16 000 (Germany only)

E-mail: info@mn-net.com

Technical support Bioanalysis

Tel.: +49 24 21 969-333

E-mail: support@mn-net.com

USA

MACHEREY-NAGEL Inc.

924 Marcon Blvd. · Suite 102 · Allentown PA, 18109 · USA

Toll-free: 888 321 6224 (MACH)

E-mail: sales-us@mn-net.com

France

MACHEREY-NAGEL SAS

1, rue Gutenberg – BP135 · 67720 Hoerdts Cedex · France

Tel.: +33 388 68 22 68

E-mail: sales-fr@mn-net.com

MACHEREY-NAGEL SAS (Société par Actions Simplifiée) au capital de 186600 €

Siret 379 859 531 00020 · RCS Strasbourg B379859531

N° intracommunautaire FR04 379 859 531

Switzerland

MACHEREY-NAGEL AG

Hirsackerstr. 7 · 4702 Oensingen · Switzerland

Tel.: +41 62 388 55 00

E-mail: sales-ch@mn-net.com