


MACHEREY-NAGEL

Dip & Read Tests

Rapid Tests



Disinfectant testing

Rapid Tests for Healthcare

- Quick & easy
- Reliable results
- Higher patient safety



MACHEREY-NAGEL

www.mn-net.com



MN Dip & Read Tests for Healthcare

Health care professionals use different types of disinfectants and sanitizers for general disinfection and for cleaning of medical devices. Prior to the disinfection of devices it proved useful to control the potency of the disinfection baths. After disinfection all devices have to be rinsed and controlled for sanitizer residues prior to the next use. Dip & Read tests allow the quick and easy detection of sanitizer residues. Thus, MN test strips increase patient safety, while making healthcare professionals' work easier and more efficient.



Product information

| Disinfectant | Product | Potency | Residual | Gradation | REF |
|--------------------|----------------------------------|---------|----------|---------------------------------------------------------------------------------------------|-------|
| Chlorine | QUANTOFIX® Chlorine 500 | ■ | | 0 · 25 · 50 · 100 · 250 · 500 mg/L Cl ₂ | 91354 |
| | QUANTOFIX® Chlorine Sensitive* | | ■ | 0 · 0.1 · 0.5 · 1 · 3 · 10 mg/L Cl ₂ | 91339 |
| Chlorine dioxide | QUANTOFIX® Chlorine dioxide 3000 | ■ | | 0 · 10 · 25 · 50 · 100 · 250 · 500 · 1000 · 1500 · 2000 · 2500 · 3000 mg/L ClO ₂ | 91357 |
| | QUANTOFIX® Chlorine dioxide 15 | | ■ | 0 · 0.25 · 0.5 · 1 · 3 · 5 · 10 · 15 mg/L ClO ₂ | 91355 |
| Glutaraldehyde | QUANTOFIX® Glutaraldehyde | ■ | | 0 · 0.5 · 1.0 · 1.5 · 2.0 · 2.5 % glutaraldehyde | 91343 |
| Hydroxyacetic Acid | pH-Fix 3.6-6.1* | | ■ | pH 3.6 · 4.1 · 4.4 · 4.7 · 5.0 · 5.3 · 5.6 · 6.1 | 92130 |
| Peracetic Acid | QUANTOFIX® Peracetic acid 2000* | ■ | | 0 · 500 · 1000 · 1500 · 2000 mg/L peracetic acid | 91342 |
| | QUANTOFIX® Peracetic acid 500* | ■ | | 0 · 50 · 100 · 200 · 300 · 400 · 500 mg/L peracetic acid | 91341 |
| | QUANTOFIX® Peracetic acid 50* | | ■ | 0 · 5 · 10 · 20 · 30 · 50 mg/L peracetic acid | 91340 |
| Peroxide | QUANTOFIX® Peroxide 1000* | ■ | | 0 · 50 · 150 · 300 · 500 · 800 · 1000 mg/L H ₂ O ₂ | 91333 |
| | QUANTOFIX® Peroxide 100* | | ■ | 0 · 1 · 3 · 10 · 30 · 100 mg/L H ₂ O ₂ | 91312 |
| | QUANTOFIX® Peroxid 25* | | ■ | 0 · 0.5 · 2 · 5 · 10 · 25 mg/L H ₂ O ₂ | 91319 |
| QUAT | QUANTOFIX® QUAT | ■ | | 0 · 10 · 25 · 50 · 100 · 250 · 500 · 1000 mg/L benzalkonium-chloride | 91337 |

* suitable for evaluation with QUANTOFIX® Relax

QUANTOFIX® Relax

The instrument is for the automatic reading of QUANTOFIX® test strips. Just place the strip on the strip adaptor and the Relax will give a precise and reliable result within seconds. The results are printed, saved internally and can be transmitted to a PC. It is ideal for QC managers and decision makers in production.

QUANTOFIX® Relax – for easy, reliable QC procedures

- highest precision
- objective results
- easy result management

REF 91346



www.mn-net.com

MACHEREY-NAGEL



MACHEREY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6-8 · 52355 Düren · Germany

DE / International:

Tel.: +49 24 21 969-0

Fax: +49 24 21 969-199

E-mail: info@mn-net.com

CH:

Tel.: +41 62 388 55 00

Fax: +41 62 388 55 05

E-mail: sales-ch@mn-net.com

FR:

Tel.: +33 388 68 22 68

Fax: +33 388 51 76 88

E-mail: sales-fr@mn-net.com

