

NANOCOLOR® gesamt-Kjeldahl-Stickstoff TKN 16
total Kjeldahl nitrogen TKN / Azote total Kjeldahl TKN / Nitrógeno total Kjeldahl NTK

Test 0-67
REF 985067

Methode / Method / Méthode / Método 345 nm
0671 1.00–32.0 mg/L TN_b
0672 1.00–16.0 mg/L NO₃-N
0673 1.00–16.0 mg/L TKN

Min	20–80 %		Max
1.00	7.20	25.8	32.0 mg/L TN _b
1.00	4.00	13.0	16.0 mg/L NO ₃ -N
1.00	4.00	13.0	16.0 mg/L TKN

1

A

2.0 mL

pH 5–9

2

2.0 mL

R1

3

4

100 °C / 60 min

5

20–25 °C

6

7

8

A

9

B

500 µL

A

10

500 µL

R2

11

3 x

12

10 min

13

14

B

M

15

C

250 µL

pH 1–13

16

250 µL

R2

17

3 x

18

10 min

19

20

C

M