


visocolor® ECO Alkalinität TA


Alkalinity TA / Alcalinité TA / Alkalinidad TA

Method(e) / Método		PF-3 Pool/ Drinking Water 450 / 590 nm	PF-12 ^{Plus} 436 / 585 nm	Advance 436 / 585 nm
5041	0.3–14.0 °d	■	■	■
5042	0.10–5.00 mmol/L H ⁺	■	■	■
5043	0.4–17.5 °e	■	■	■
5044	0.5–25.0 °f	■	■	■
5045	5–250 mg/L CaCO ₃	■	■	■




Null mit Probe
Blank with sample
Blanc avec échantillon
En blanco con muestra


Geeignet für Meerwasser
Suitable for seawater
Convient à l'eau de mer
Adecuado para el agua de mar




Küvetten spülen
Rinse test tube
Rincer les éprouvettes
Lavar los tubos




Probe
Sample
Echantillon
Muestra




200 µL TA-1



1 NANOFIX TA-2




Schütteln
Shake
Agiter
Agitar



Säubern
Clean
Nettoyer
Limpiar



Warten
Wait
Attendre
Espere



Messung
Measurement
Mesure
Medición

Umrechnungstabelle / Conversion table / Tableau de conversion / Tabla de conversión

°d	°e	°f	mg/L CaCO ₃	mmol/L H ⁺	gpg
1	1.3	1.8	18	0.36	1
2	2.5	3.6	36	0.72	2
3	3.8	5.4	54	1.08	3
4	5.0	7.1	71	1.42	4
5	6.3	8.9	89	1.78	5
6	7.5	10.7	107	2.14	6
7	8.8	12.5	125	2.50	7
8	10.0	14.3	143	2.86	8
9	11.3	16.1	161	3.22	9
10	12.5	17.8	178	3.56	10