


visicolor[®] ECO Cyanursäure


Cyanuric acid / Acide cyanurique / Ácido cianúrico

| Method(e) / Método | PF-3 Pool/ Drinking Water 450 nm | PF-12 540 nm | PF-12 ^{Plus} 540 nm | Advance 540 nm |
|----------------------|--|-----------------|---------------------------------|-------------------|
| 5231 10–100 mg/L Cya | ■ | ■ | ■ | ■ |




Null mit Probe
Blank with sample
Blanc avec échantillon
En blanco con muestra


Geeignet für Meerwasser
Suitable for seawater
Convient à l'eau de mer
Adecuado para el agua de mar




Küvetten spülen
Rinse test tube
Rincer les éprouvettes
Lavar los tubos




Probe
Sample
Echantillon
Muestra



1/2 Cya-1




Schütteln
Shake
Agiter
Agitar



Säubern
Clean
Nettoyer
Limpiar



Warten
Wait
Attendre
Espere



Messung
Measurement
Mesure
Medición