


visicolor® ECO Kupfer


Copper / Cuivre / Cobre

Method(e) / Método


5371 0.1–5.0 mg/L Cu²⁺PF-3 Fish
530 nmPF-12
585 nmPF-12^{Plus}
585 nmAdvance
600 nm

Null mit Probe
Blank with sample
Blanc avec échantillon
En blanco con muestra


Geeignet für Meerwasser
Suitable for seawater
Convient à l'eau de mer
Adecuado para el agua de mar




Küvetten spülen
Rinse test tube
Rincer les éprouvettes
Lavar los tubos




Probe
Sample
Echantillon
Muestra




5 ♂ Cu-1



Schütteln
Shake
Agiter
Agitar




5 ♂ Cu-2




Schütteln
Shake
Agiter
Agitar



Säubern
Clean
Nettoyer
Limpiar



Warten
Wait
Attendre
Espere



Messung
Measurement
Mesure
Medición